




Tecnológico
de Monterrey



Rheometry




Questions and Reflections

- Besides of the capillary test,
are there some other ways to measure the viscosity?
- What other moduli exist for solid and liquids?
- Are those moduli dependent from each other?
- If they are dependent, how can a given modulus be calculated from another one?

Melt flow indexer

- One of the most used equipment in the polymer industry
- A quality control tool
- Measures “flowability”
 - Melt flow rate: # grams/10 minutes
- According to Standards:
 - ISO 1133 and ASTM D1238



Brookfield Viscometer

- For each “spindle” velocity, gives the modulus viscosity, based on the torque.
- The viscosity is typically reported in:
 - Centipoise or mPa.s




Other rheological tests:

- **Constant shear rate** (viscous modulus)
- **Constant stress** (creep and recovery compliance)
- **Constant strain** (relaxation modulus)



VISCOUS MODULUS

Viscous modulus: the material's response to a constant shear rate



$$\text{MODULUS} = \frac{\text{STRESS}}{\text{SHEAR RATE}}$$

Viscous modulus

$$\eta = \sigma / \dot{\gamma}$$

$$\dot{\gamma} = dv_{\theta} / dz$$


σ in Pa

and


shear rate in 1/s

η in Pa.s


Viscosity versus time for 4 different HDPE resins when a sudden shear rate of 0.01 1/s is imposed in a parallel plate rheometer



Viscosity versus time (Capillary rheometer)




Shear viscosity versus shear rate for a polypropylene resin



Shear viscosity vs. shear rate


(steady state shear viscosity vs. constant shear rate)
for a set of reactor grade polypropylene resins





CONSTANT STRESS RHEOMETRY

Constant stress rheometry: creep and recovery compliance



$$\text{MODULUS} = \text{STRAIN} / \text{STRESS}$$

Compliance modulus

$$J = \gamma / \sigma$$


$$\gamma = (L - L_0) / L_0$$

σ in Pa

then

J in 1/Pa


Constant stress rheology:



But, once the stress ceases the material...

Recoils and the extent of recoil with respect to the extent of deformation under constant stress is called the the **recovery compliance**


Constant stress rheology: creep and recovery compliance




Where:

- J** is the compliance: Strain / Stress
- J_r (0)** is the compliance at the time at which the stress ceases (t_c)
- J_r** is the compliance at any time after the stress ceases.
- J_r(0)-J_r** is the recoverable compliance at a given time ($t - t_c$)
- J_e⁰** is the steady state recoverable compliance
 $J_r(0) - J_r(t_{\text{infinite}})$ (at $t \gg t_c$).
- J_g** instantaneous or glassy compliance.

Creep compliance vs. time (strain under constant stress vs time) for a set of reactor grade polypropylene resins



Recovery compliance vs. time (strain under constant stress vs time) for a set of reactor grade polypropylene resins



Creep and recovery compliance data for the RGPP resins


applied constant stress: 300 Pa

resin	time	strain	$J_{\max} \times 10^{-5}$ at 300 s	shear rate	eta	$J_e^0 \times 10^{-4}$ at 300 s
par0	316	4.05	1351	1.24E-02	24200	4.95
pbr0	316	16.98	5660	5.34E-02	5620	3.89
pcr0	316	19.95	6650	6.28E-02	4780	3.86
pdr0	316	35.70	11900	1.11E-01	2710	3.51
pcont	316	30.63	10210	9.55E-02	3140	4.41



RELAXATION MODULUS

Relaxation modulus: the material's response to a sudden strain



$$\text{MODULUS} = \text{STRESS} / \text{STRAIN}$$

Relaxation modulus

$$G = \sigma / \gamma$$

$$\gamma = (L - L_0) / L_0$$


σ in Pa

then

G in Pa


Typical Relaxation Moduli

Relaxation moduli at different % strain for a HDPE resin at
 $T=210^{\circ}\text{C}$



Typical Relaxation Moduli

Relaxation Moduli for HDPE resins after a sudden 100% strain at $T=210^{\circ}\text{C}$





Tecnológico de Monterrey

Derechos Reservados 2015 Tecnológico de Monterrey

Prohibida la reproducción total o parcial de esta obra sin expresa autorización del Tecnológico de Monterrey.