

Characterization techniques applied to Carbon Nanowires

Team 3

Angel Manuel Villalba Rodríguez <A00828035@itesm.mx>

Constanza Alvarez López <A00829469@itesm.mx>

Antonio Osamu Katagiri Tanaka <A01212611@itesm.mx>

Thursday, 04 June 2020




Agenda

- Project Overview, goals & justification : <2min
- Characterization Techniques : ~10min
 - Rheometry
 - Scanning Electron Microscopy (SEM)
 - Raman Spectroscopy
 - 4-Wire Kelvin
- Conclusions : <1min
- Q&A : <4min



Justification:




J. Jang, J. Bae, Carbon nanofiber/polypyrrole nanocable as toxic gas sensor, Sensors Actuators B Chem. 122 (2007) 7–13. <https://doi.org/10.1016/j.snb.2006.05.002>.

X. Shi, W. Zhou, D. Ma, Q. Ma, D. Bridges, Y. Ma, A. Hu, Electrospinning of Nanofibers and Their Applications for Energy Devices, J. Nanomater. 2015 (2015) 1–20. <https://doi.org/10.1155/2015/140716>.



C. Monterrey, B. Cárdenas Benítez, INSTITUTO TECNOLÓGICO Y DE ESTUDIOS SUPERIORES DE MONTERREY Advanced Manufacturing Techniques for the Fabrication and Surface Modification of Carbon Nanowires, (2017) 160.

Introduction: Project overview

To design polymer solutions that can be **electrospun** by NFES, **photopolymerized**, and then **pyrolyzed** into conductive carbon nano-wires.



Charc. Techniques: Rheology (instrumentation)



Measured Parameter	Calculated Parameter
Torque	Stress
Angular displacement	Strain
Angular velocity	Shear rate

$$\text{Stress} = \text{Torque}[\text{N} \cdot \text{m}] \cdot \text{Stress constant}$$

$$\text{Strain} = \text{Angular displacement}[\text{rad}] \cdot \text{Strain constant}$$

$$\text{Shear rate} = \text{Angular velocity} \left[\frac{\text{rad}}{\text{s}} \right] \cdot \text{Strain constant}$$

$$\text{Viscosity} = \frac{\text{Stress}}{\text{Shear rate}}$$

$$\text{Young's Modulus} = \frac{\text{Stress}}{\text{Strain}}$$



Charc. Techniques: Rheology (process parameters)

Important rheometer parameters:

- Torque Range
- Angular Resolution
- Angular Velocity Range
- Frequency Range


Geometry options:

Concentric
Cylinders

Cone and
Plate

Parallel
Plate


Torsion
Rectangular



Very Low Viscosity → Solids




Charc. Techniques: Rheometry (Example results)



Charc. Techniques: 4-Wire Kelvin


Typical method:



$$\text{Ohmmeter}_{\text{indication}} = R_{\text{wire1}} + R_{\text{subject}} + R_{\text{wire2}}$$

If the sample has a very low resistance, the measurement **error introduced by wire resistance** will be substantial.

Kelvin method:



$$R_{\text{subject}} = \frac{\text{Voltmeter}_{\text{indication}}}{\text{Ammeter}_{\text{indication}}}$$

Electrical conductivity:

$$\sigma = \frac{L}{R_{\text{subject}} \cdot A} \left[\frac{S}{m} \right]$$




Charc. Techniques: Scanning Electron Microscopy (SEM)


A beam of electrons scans across the material surface, producing high-energy back-scattered electrons and low-energy secondary electrons. The back-scattered electrons provides a physical image of the sample surface.

Sample preparation:
After assesing conductivity of carbon nanowires, **metal coating** is only applied when they **don't present conductivity**

Fiber diameter **before** pyrolysis



Fiber diameter **after** pyrolysis



Charc. Techniques: Raman Spectroscopy

Collect knowledge about **molecular information** and **crystalline structures**. A laser irradiate the sample with monochromatic radiation, with excitations in the UV, visible and near-IR spectral region (785-1064nm). The irradiation generates a Raman spectrum of scattered light.


Sample preparation:

The pyrolyzed sample is placed in a **quartz or glass vessel** and then analyzed, no more sample preparation is needed for Raman.



Charc. Techniques: Raman Spectroscopy

Expectation is to find two bands (marked D and G in the figure), centered about 1360 and 1594 cm^{-1} to confirm the presence of graphitic carbon in the nanofibers.



Conclusions:

- The presented techniques provide enough data to satisfy the project goal.

Technique	Data obtained
Rheology	Zero-shear viscosity
4-Wire Kelvin	Electrical resistivity
SEM	Fiber morphology
Raman	Presence of graphitic carbon

- The Rheology ensures electrospunable solutions
- TGA may provide useful information, but is omitted due to time and availability constraints.



Characterization techniques applied to Carbon Nanowires

Team 3

Angel Manuel Villalba Rodríguez <A00828035@itesm.mx>

Constanza Alvarez López <A00829469@itesm.mx>

Antonio Osamu Katagiri Tanaka <A01212611@itesm.mx>

Thursday, 04 June 2020

